



Summary of the most commonly used test methods

Further Total Joint Replacement Implants

Total ankle-joint replacement



ISO 22622 - Implants for surgery — Wear of total ankle-joint prostheses — Loading and displacement parameters for wear-testing machines with load or displacement control and corresponding environmental conditions for test




ASTM F2665 - Standard specification for total ankle replacement prosthesis



ISO 21534 - Non-active surgical implants — Joint replacement implants — Particular requirements

Total elbow replacement



ASTM F2887 - Standard specification for total elbow prostheses



ISO 21534 - Non-active surgical implants — Joint replacement implants — Particular requirements

Small joint replacement



ISO 21534 - Non-active surgical implants — Joint replacement implants — Particular requirements

For further information please contact:

SpineServ GmbH & Co. KG, Soefflinger Strasse 100, 89077 Ulm, Germany

Phone: +49 (0)731 175 6788

experts@spineserv.de

www.spineserv.com