


CATHETERS – Test Methods


Enteral feeding systems

ASTM F2528




Intravascular introducers, dilators, guidewires

ISO 11070




Catheter systems

DIN EN 13868



ISO 20697



- ISO 20695
- ISO 20696
- ISO 20698

Intravascular catheters

ISO 10555 series of standards

- ISO 10555-1
 - ISO 10555-2
 - ISO 10555-3
 - ISO 10555-4
 - ISO 10555-5
 - ISO 10555-6
- » General requirements
 - » Angiographic catheters
 - » Central venous catheters
 - » Balloon dilatation catheters
 - » Over-needle peripheral catheters
 - » Subcutaneous implanted ports

